



# Combined Electrocardiography- and Respiratory-Triggered CT of the Lung to Reduce Respiratory Misregistration Artifacts between Imaging Slabs in Free-Breathing Children: Initial Experience

Hyun Woo Goo, MD, PhD<sup>1</sup>, Thomas Allmendinger, PhD<sup>2</sup>

<sup>1</sup>Department of Radiology and Research Institute of Radiology, University of Ulsan College of Medicine, Asan Medical Center, Seoul 05505, Korea; <sup>2</sup>Siemens Healthcare, GmbH, Computed Tomography Division, Forchheim 91301, Germany

**Objective:** Cardiac and respiratory motion artifacts degrade the image quality of lung CT in free-breathing children. The aim of this study was to evaluate the effect of combined electrocardiography (ECG) and respiratory triggering on respiratory misregistration artifacts on lung CT in free-breathing children.

**Materials and Methods:** In total, 15 children (median age 19 months, range 6 months–8 years; 7 boys), who underwent free-breathing ECG-triggered lung CT with and without respiratory-triggering were included. A pressure-sensing belt of a respiratory gating system was used to obtain the respiratory signal. The degree of respiratory misregistration artifacts between imaging slabs was graded on a 4-point scale (1, excellent image quality) on coronal and sagittal images and compared between ECG-triggered lung CT studies with and without respiratory triggering. A  $p$  value  $< 0.05$  was considered significant.

**Results:** Lung CT with combined ECG and respiratory triggering showed significantly less respiratory misregistration artifacts than lung CT with ECG triggering only ( $1.1 \pm 0.4$  vs.  $2.2 \pm 1.0$ ,  $p = 0.003$ ).

**Conclusion:** Additional respiratory-triggering reduces respiratory misregistration artifacts on ECG-triggered lung CT in free-breathing children.

**Keywords:** Lung CT; ECG triggering; Respiratory triggering; Children; Motion artifacts

## INTRODUCTION

It is desirable to minimize motion artifacts from cardiac pulsation and respiration on lung CT in free-breathing children. To achieve this, faster gantry rotation time (1), partial reconstruction (2), electrocardiography (ECG) triggering (1, 3, 4), respiratory triggering (5–8), controlled

ventilation (9–11), and general anesthesia with intubation (12) have been used. The use of faster gantry rotation time and partial reconstruction is insufficient to suppress all motion artifacts on free-breathing lung CT (13). The respiratory-controlled techniques require intense patient cooperation or handling by skilled experts; hence, they are difficult to apply in free-breathing children (5, 7, 8–12). In contrast, transthoracic impedance plethysmography-based respiratory triggering used in a study (6) may be applicable in free-breathing children. Retrospective respiratory gating has been used for radiotherapy (14), positron emission tomography (PET)/CT (15), and animal lung CT (16, 17). However, a high radiation dose of retrospective respiratory-gated lung CT hinders its clinical use in children.

Combined ECG and respiratory triggering or gating is used to suppress motion artifacts during free breathing in lung MRI (18), or cardiac PET (19). For free-breathing lung CT, a combined ECG and respiratory triggering system may also

Received December 4, 2016; accepted after revision March 12, 2017.

**Corresponding author:** Hyun Woo Goo, MD, PhD, Department of Radiology and Research Institute of Radiology, University of Ulsan College of Medicine, Asan Medical Center, 88 Olympic-ro, 43-gil, Songpa-gu, Seoul 05505, Korea.

• Tel: (822) 3010-4388 • Fax: (822) 476-0090  
• E-mail: hwgoo@amc.seoul.kr

This is an Open Access article distributed under the terms of the Creative Commons Attribution Non-Commercial License (<http://creativecommons.org/licenses/by-nc/4.0>) which permits unrestricted non-commercial use, distribution, and reproduction in any medium, provided the original work is properly cited.

be used to minimize respiratory misregistration artifacts between adjacent imaging slabs, which is frequently seen with ECG-triggering alone; however, to our best knowledge, its technical feasibility has not been evaluated previously. Prospective ECG-triggered sequential scanning is already widely available in clinical CT scanners (20, 21) and additional prospective respiratory triggering, hence, needs to be incorporated for a combined triggering system. Therefore, the aim of this study was to describe our initial experience with combined ECG and respiratory triggering to reduce respiratory misregistration artifacts on lung CT in free-breathing children.

## MATERIALS AND METHODS

### Study Population


Our Institutional Review Board approved the retrospective study; and all CT examinations were clinically indicated. In our institution, prospectively ECG-triggered sequential scanning is used for unenhanced lung CT in free-breathing sedated children to eliminate cardiac motion artifacts and reduce respiratory motion artifacts on axial images. However, respiratory misregistration artifacts persist on coronal and sagittal reformatted images between adjacent sequential imaging slabs. In addition to axial images, coronal and sagittal reformatted images are included in routine image interpretation; hence, there is a need to reduce the motion artifacts on reformatted images. In our institution, lung CT protocol was gradually changed from prospectively ECG-triggered sequential scanning to combined prospectively ECG- and respiratory-triggered sequential scanning to reduce the respiratory misregistration artifacts on reformatted images. Combined prospectively ECG- and respiratory-triggered scan mode was installed in our institution in December 2008. Depending on the patients' condition and examination schedule, prospectively ECG-triggered sequential or spiral scan was used for unenhanced chest CT until March 2011. After March 2011, combined prospectively ECG- and respiratory-triggered scan was used as the primary scan mode for unenhanced chest CT in free-breathing children. In this study, we included pediatric patients who consecutively underwent free-breathing unenhanced lung CT using these two different CT scan protocols. Patients older than 15 years of age were excluded.

In 15 children (median age 19 months, range 6 months–8 years; 7 boys), prospectively ECG-triggered sequential


lung CT using a first-generation dual-source scanner (SOMATOM Definition, Siemens Healthcare, Forchheim, Germany) was obtained during free breathing with and without additional respiratory triggering on two different days (median interval 3 months, range 12 days–1 year) between December 2008 and March 2011. Lung CT with ECG triggering was performed first in 10 patients and later in five patients. Disease categories of the patient population comprised congenital pulmonary airway malformation (n = 7), pneumonia (n = 2), Wilms tumor (n = 2), interstitial pneumonitis (n = 1), Langerhans cell histiocytosis (n = 1), repaired total anomalous pulmonary venous connection (n = 1), and pectus excavatum (n = 1). All children were sedated with oral chloral hydrate and intravenous administration of midazolam and/or ketamine was added if necessary. All children were younger than 5 years of age except for an 8-year-old boy who needed sedation for lung CT due to altered mentality caused by hemorrhagic brain lesions. In all children, general anesthesia, intubation, or controlled ventilation was not used.

### CT Imaging

Dual-source prospective ECG-triggered sequential CT scan was performed at a gantry rotation time of 0.33 second (temporal resolution of 83 ms), 2 x 32 x 0.6 mm detector collimation with a z-flying focal spot, 0.75-mm slice width, 0.4-mm reconstruction interval, and a cycle time of 1.36 seconds for one acquisition and the subsequent table feed resulting in data acquisitions every 2 or 3 R-R intervals depending on heart rates. Sequential scanning was acquired at the end-systole of the R-R interval by using an absolute trigger delay (ms). For combined ECG and respiratory triggering, a pressure-sensing belt of a respiratory gating system (AZ-733V, Anzai Medical Co., Tokyo, Japan) was placed, usually in the upper abdomen, to obtain the respiratory signal (Fig. 1). Then, respiratory triggering was produced at either inspiration (n = 5) or expiration peak (n = 10) of the respiratory cycle determined by a pediatric radiologist based on the shape of respiratory curve in each child. Sequential scanning was acquired only when ECG triggering was synchronous with respiratory triggering through a prototype box that was not commercially available (Fig. 1). Average sequential scan rates per minute during free-breathing lung CT were recorded and compared between the two scan modes. Body size-adapted volume CT dose index derived from 16 cm-sized phantom of lung CT was determined individually by the cross-sectional area



**Fig. 1. Schematic drawing of combined ECG and respiratory triggering.** From ECG leads, R wave (dot) is detected and end-systolic phase (arrows) is determined by trigger delay. Pressure-sensing belt is placed in upper abdomen to continuously monitor respiratory signal; in this case, respiratory triggering (dot) is produced at inspiration peak. When they occur concurrently, these two triggering signals can be passed through gate box to CT system. Combined ECG and respiratory triggering at R wave on ECG (dot) is displayed on CT console. From combined ECG and respiratory triggering, end-systolic sequential scanning is acquired every three R waves in this case. ECG = electrocardiography



**Fig. 2. Measurement of respiratory misregistration artifacts between imaging slabs on sagittal formatted lung CT image.** Respiratory misregistration occurred at three levels and degree of misregistration was measured in length (brackets) at six locations. Among them, maximal length (Max, orange yellow bracket) was used for semi-quantitative grading (grade 1, < 2 mm; grade 2, 2–5 mm; grade 3, 5–10 mm; and grade 4, > 10 mm).

and the mean body density of each patient (22). To reduce image noise of unenhanced lung CT images, high tube voltage was used, e.g., 100 kVp in 14 patients, 120 kVp in one patient. Multiple sequential scans were obtained to cover the whole lung. The number of sequential scan slabs per lung CT was compared between ECG-triggered lung CT and combined ECG- and respiratory-triggered lung CT. CT dose estimates were calculated from dose length products derived from 16 cm-sized CT dose index phantom and age-, sex-, and tube voltage-specific conversion factors (23). The calculated effective dose was compared between ECG-triggered lung CT and combined ECG- and respiratory-triggered lung CT. For image interpretation, 3- to 5-mm thick axial, coronal, and sagittal CT images were reconstructed depending on scan length.

**Image Evaluation**

The degree of respiratory misregistration artifacts in the entire thorax including the chest wall, the lung outline, and the diaphragm was measured by using calipers on coronal

and sagittal lung CT images (Fig. 2) and subsequently, semi-quantitatively graded on a 4-point scale (grade 1, maximal distance of misalignment between adjacent imaging slabs due to respiratory misregistration < 2 mm; grade 2, 2–5 mm; grade 3, 5–10 mm; and grade 4, > 10 mm).

**Statistical Analysis**


For statistical analysis, statistical software (SPSS version

20.0, IBM, Chicago, IL, USA) was used. Variables were expressed as the mean with standard deviation. Respiratory misregistration artifact grades were compared between ECG-triggered lung CT studies with and without respiratory triggering by using the Wilcoxon signed rank test. Continuous variables were compared between the two lung CT studies by using the paired *t* test. A *p* value < 0.05 was considered as statistically significant.

**Table 1. Comparisons of Variables between ECG-Triggered Lung CT with and without Respiratory Triggering**

	ECG Triggering Only	Combined ECG and Respiratory Triggering	<i>P</i>
Respiratory misregistration artifact grade	2.2 ± 1.0	1.1 ± 0.4	0.003
Scan rates per minute	50.2 ± 6.8	35.3 ± 18.9	0.01
Imaging slab number	8.3 ± 1.9	8.1 ± 1.8	0.6
Effective dose (mSv)	1.0 ± 0.3	1.1 ± 0.3	0.2

ECG = electrocardiography



**Fig. 3. 5-year-old boy with Langerhans cell histiocytosis.**


**A, B.** Free-breathing coronal (A) and sagittal (B) ECG-triggered lung CT images acquired at end-systole show grade 4 respiratory misregistration artifacts. **C, D.** Respiratory misregistration artifacts grade was grade 1 on free-breathing coronal (C) and sagittal (D) lung CT images with combined ECG and respiratory triggering at end-systole and inspiration peak.

## RESULTS

Combined triggering was successfully applied to all children without technical failures. The results of comparisons between ECG-triggered lung CT studies with and without respiratory triggering are described in Table 1. Respiratory misregistration artifact grades on lung CT were significantly less with combined ECG and respiratory triggering ( $1.1 \pm 0.4$ ) than with ECG-triggering only ( $2.2 \pm 1.0$ ,  $p = 0.003$ ) (Figs. 3, 4). Average sequential scan rates per minute were significantly higher for ECG-triggered scan than for combined ECG- and respiratory-triggered scan ( $50.2 \pm 6.8$  vs.  $35.3 \pm 18.9$ ,  $p = 0.01$ ). The number of sequential scan slabs per lung CT was not significantly different between ECG triggering only ( $8.3 \pm 1.9$ ) and combined ECG and respiratory triggering ( $8.1 \pm 1.8$ ). The effective CT dose values showed no significant differences between ECG-triggered lung CT studies with and without respiratory triggering ( $1.0 \pm 0.3$  mSv vs.  $1.1 \pm 0.3$  mSv,  $p = 0.2$ ).

## DISCUSSION

Image quality of lung CT in free-breathing children has improved mainly due to more rapid acquisition of CT data. However, residual motion-related artifacts need to be reduced to achieve comparable image quality on free-breathing lung CT as that on breath-hold lung CT. ECG-triggered sequential scanning is often used to reduce cardiac pulsation artifacts on lung CT, but respiratory misregistration artifacts are particularly pronounced along the z-axis in free-breathing children. In this study, we demonstrated the technical feasibility of combined ECG- and respiratory-triggered lung CT to minimize not only cardiac pulsation artifacts, but also respiratory misregistration artifacts in free-breathing children without any interventions to patients' airways, such as general anesthesia with intubation and controlled ventilation. In the 1990s, the first respiratory-triggered lung CT using a spirometer was reported in adult patients (5, 7). However, it has not been widely used since the spirometer is not patient-friendly. Recently, Lell et al. (24) reported that free-breathing respiratory-triggered high-pitch dual-source chest spiral CT, using the same respiratory gating system as in this study, could minimize motion artifacts in adult patients, as compared with conventional inspiratory breath-hold spiral chest CT. However, the sole effect of respiratory triggering on reduction of motion artifacts could not be



**Fig. 4.** 2-year-old girl with multifocal congenital pulmonary airway malformation.

Free-breathing axial (A), coronal (B), and sagittal (C) lung CT images with combined ECG and respiratory triggering at end-systole and expiration peak demonstrate excellent normal anatomic and pathologic details without any motion artifacts.

determined, because high-pitch dual-source spiral scan itself can considerably reduce motion artifacts due to its very fast scan time. To our best knowledge, combined ECG- and respiratory-triggered lung CT has not been previously evaluated in free-breathing patients.

In the combined ECG- and respiratory-triggered lung CT, sequential scan rates were approximately 1.4-times slower than those in ECG-triggered lung CT, which prolonged the total scan time of combined-triggered scan. However, this prolongation may not affect the image quality of lung CT without intravenous contrast agent. In contrast, a modified intravenous contrast injection protocol to lengthen peak enhancement period is necessary for contrast-enhanced chest or cardiac CT with combined ECG and respiratory triggering in free-breathing patients. This can be achieved by using multiphase decremental injection rates of contrast agent (25, 26).

In addition to improved anatomic details, the combined ECG and respiratory triggering used in this study has potential for functional and quantitative lung evaluations, such as air trapping evaluation and lung volumetry, in free-breathing children because motion-free lung CT image data can be obtained at both end-inspiration and end-expiration with three-dimensional non-rigid image registration (27, 28). The combined ECG and respiratory triggering can also be used for uncooperative, free-breathing adults with altered mental state or respiratory difficulty. The respiratory gating system used in this study has a technical limitation, since respiratory triggering is not possible when the respiratory rate is higher than 50/min or breathing is substantially irregular.

This study has several limitations. First, the study population was small. However, the results of paired comparisons to evaluate the technical feasibility in this retrospective study may not be different even by increasing the number of patients included. Second, underlying diseases of the study population were heterogeneous. However, this was a technical feasibility study focusing on motion artifacts rather than clinical diseases. Third, there was a median 3-months-time interval between the two lung CT studies in the same patient. This was also attributed to the retrospective nature using the clinically indicated CT studies.

In conclusion, combined ECG and respiratory triggering can minimize respiratory misregistration artifacts between imaging slabs on lung CT in free-breathing children. This new CT imaging technique may eliminate the need for

more invasive procedures, such as general anesthesia with intubation and controlled ventilation, to produce a motion-free effect on free-breathing lung CT.

#### Conflict of Interest

Thomas Allmendinger is a Siemens Healthcare employee.

#### REFERENCES

1. Montaudon M, Berger P, Blachère H, De Boucaud L, Latrabe V, Laurent F. Thin-section CT of the lung: influence of 0.5-s gantry rotation and ECG triggering on image quality. *Eur Radiol* 2001;11:1681-1687
2. Ha HI, Goo HW, Seo JB, Song JW, Lee JS. Effects of high-resolution CT of the lung using partial versus full reconstruction on motion artifacts and image noise. *AJR Am J Roentgenol* 2006;187:618-622
3. Schoepf UJ, Becker CR, Bruening RD, Helmberger T, Staebler A, Leimeister P, et al. Electrocardiographically gated thin-section CT of the lung. *Radiology* 1999;212:649-654
4. Nishiura M, Johkoh T, Yamamoto S, Honda O, Kozuka T, Koyama M, et al. Electrocardiography-triggered high-resolution CT for reducing cardiac motion artifact: evaluation of the extent of ground-glass attenuation in patients with idiopathic pulmonary fibrosis. *Radiat Med* 2007;25:523-528
5. Kalender WA, Rienmüller R, Seissler W, Behr J, Welke M, Fichte H. Measurement of pulmonary parenchymal attenuation: use of spirometric gating with quantitative CT. *Radiology* 1990;175:265-268
6. Mori M, Murata K, Takahashi M, Shimoyama K, Ota T, Morita R, et al. Accurate contiguous sections without breath-holding on chest CT: value of respiratory gating and ultrafast CT. *AJR Am J Roentgenol* 1994;162:1057-1062
7. Robinson TE, Leung AN, Moss RB, Blankenberg FG, al-Dabbagh H, Northway WH. Standardized high-resolution CT of the lung using a spirometer-triggered electron beam CT scanner. *AJR Am J Roentgenol* 1999;172:1636-1638
8. Zaporozhan J, Ley S, Unterhinninghofen R, Saito Y, Fabel-Schulte M, Keller S, et al. Free-breathing three-dimensional computed tomography of the lung using prospective respiratory gating: charge-coupled device camera and laser sensor device in an animal experiment. *Invest Radiol* 2006;41:468-475
9. Long FR, Castile RG, Brody AS, Hogan MJ, Flucke RL, Filbrun DA, et al. Lungs in infants and young children: improved thin-section CT with a noninvasive controlled-ventilation technique--initial experience. *Radiology* 1999;212:588-593
10. Long FR, Castile RG. Technique and clinical applications of full-inflation and end-exhalation controlled-ventilation chest CT in infants and young children. *Pediatr Radiol* 2001;31:413-422
11. Mueller KS, Long FR, Flucke RL, Castile RG. Volume-monitored chest CT: a simplified method for obtaining motion-free images near full inspiratory and end expiratory lung volumes. *Pediatr Radiol* 2010;40:1663-1669

12. Lee EY, Mason KP, Zurakowski D, Waltz DA, Ralph A, Riaz F, et al. MDCT assessment of tracheomalacia in symptomatic infants with mediastinal aortic vascular anomalies: preliminary technical experience. *Pediatr Radiol* 2008;38:82-88
13. Ritchie CJ, Godwin JD, Crawford CR, Stanford W, Anno H, Kim Y. Minimum scan speeds for suppression of motion artifacts in CT. *Radiology* 1992;185:37-42
14. Thiyagarajan R, Sinha SN, Ravichandran R, Samuvel K, Yadav G, Sigamani AK, et al. Respiratory gated radiotherapy-pretreatment patient specific quality assurance. *J Med Phys* 2016;41:65-70
15. Werner MK, Parker JA, Kolodny GM, English JR, Palmer MR. Respiratory gating enhances imaging of pulmonary nodules and measurement of tracer uptake in FDG PET/CT. *AJR Am J Roentgenol* 2009;193:1640-1645
16. Ley S, Ley-Zaporozhan J, Unterhinninghofen R, Saito Y, Fabel-Schulte M, Weinheimer O, et al. Investigation of retrospective respiratory gating techniques for acquisition of thin-slice 4D-multidetector-computed tomography (MDCT) of the lung: feasibility study in a large animal model. *Exp Lung Res* 2006;32:395-412
17. Behzadi C, Groth M, Henes FO, Schwarz D, Deibele A, Begemann PG, et al. Intraindividual comparison of image quality using retrospective and prospective respiratory gating for the acquisition of thin sliced four dimensional multidetector CT of the thorax in a porcine model. *Exp Lung Res* 2015;41:489-498
18. Oechsner M, Pracht ED, Staeb D, Arnold JF, Köstler H, Hahn D, et al. Lung imaging under free-breathing conditions. *Magn Reson Med* 2009;61:723-727
19. Teräs M, Kokki T, Durand-Schaefer N, Noponen T, Pietilä M, Kiss J, et al. Dual-gated cardiac PET-clinical feasibility study. *Eur J Nucl Med Mol Imaging* 2010;37:505-516
20. Goo HW, Yang DH. Coronary artery visibility in free-breathing young children with congenital heart disease on cardiac 64-slice CT: dual-source ECG-triggered sequential scan vs. single-source non-ECG-synchronized spiral scan. *Pediatr Radiol* 2010;40:1670-1680
21. Goo HW. Current trends in cardiac CT in children. *Acta Radiol* 2013;54:1055-1062
22. Goo HW. Individualized volume CT dose index determined by cross-sectional area and mean density of the body to achieve uniform image noise of contrast-enhanced pediatric chest CT obtained at variable kV levels and with combined tube current modulation. *Pediatr Radiol* 2011;41:839-847
23. Goo HW. CT radiation dose optimization and estimation: an update for radiologists. *Korean J Radiol* 2012;13:1-11
24. Lell MM, Scharf M, Eller A, Wuest W, Allmendinger T, Fuchs F, et al. Feasibility of respiratory-gated high-pitch spiral CT: free-breathing inspiratory image quality. *Acad Radiol* 2016;23:406-412
25. Bae KT. Peak contrast enhancement in CT and MR angiography: when does it occur and why? Pharmacokinetic study in a porcine model. *Radiology* 2003;227:809-816
26. Goo HW. State-of-the-art CT imaging techniques for congenital heart disease. *Korean J Radiol* 2010;11:4-18
27. Goo HW. State-of-the-art pediatric chest imaging. *Pediatr Radiol* 2013;43:261
28. Goo HW. Advanced functional thoracic imaging in children: from basic concepts to clinical applications. *Pediatr Radiol* 2013;43:262-268